

## GLUCOSE - 6 - PHOSPHATE DEHYDROGENASE

This reagent kit is for qualitative estimation of **GLUCOSE-6 PHOSPHATE DEHYDROGENASE** activity in erythrocytes.

### PRINCIPLE:

The substrate Glucose-6-Phosphate is oxidized to 6-phosphogluconate with the reduction NADP to NADPH. Reaction is catalysed by Glucose-6- Phosphate dehydrogenase that is present in erythrocytes. Reduction of NADP to NADPH decolourises the blue coloured indophenol dye (DCPIP) to a brick red colour which is contributed by the hemolysate. This reagent uses a potent inhibitor to 6-Phosphogluconate dehydrogenase thus preventing its interference, in the estimation of G6PD



### CLINICAL SIGNIFICANCE :

Estimating G6PDH activity is a first important step in screening hemolytic episodes and anemia. Hereditary abnormalities can be the cause of enzyme deficiencies.

Subject deficient in this enzyme are prone to hemolytic episodes and hemolytic anemias when treated with certain anti-malaria or sulfa drugs.

### REAGENTS:

All the reagents are to be stored at 2 -8° c

	No. of bottles	
	8x0.5ml	15x0.5ml
Reagent 1 (Enzyme/coenzyme)	8	15
Reagent 2 (Buffer)	1	1
Reagent 3 (Lysing Reagent)	1	1
Reagent 4 (Inert oil)	1	1

### SPECIMEN COLLECTION AND STORAGE :

- \* Fresh whole blood is necessary.
- \* Sample should be collected preferably in EDTA as an Anticoagulant.
- \* The hemoglobin content of the sample should be Estimated.

### PRECAUTION:

- \* Estrom G-6-PDH reagent is for In Vitro diagnostic use Only.

### REAGENT RECONSTITUTION:

(A) 8x0.5 ml & 15x0.5 ml: Reconstitute 1 vial of Reagent 1(Enzyme/Chromogen) using 0.5 ml of Reagent 2 (Buffer). The reagent should be reconstituted just before use and should be protected from light.

### REAGENT STORAGE AND STABILITY:

All the reagents are stable up to expiry date stated on the label.

### PROCEDURE:

Allow the sample and reagent to attain room temperature prior to use. Determine the hemoglobin content of the blood Sample

Dispense into test tubes	Test
Reagent 3 (Precooled Lysing Reagents)	1.0 ml
Sample (Fresh Whole Blood)	20 µl

Mix well and keep at 2-8° c for 5 minutes.

**Note:** if the hemoglobin content of the sample is significantly less/more than 15gms/dl, the sample volume may be adjusted as follows.

Hemoglobin Concentration	Sample Volume
13 gms/dl	23 µl
14 gms/dl	21 µl
16 gms/dl	19 µl

### PROCEDURE:

1. The hemolysate prepared as stated above should be transferred completely to the freshly prepared working reagents Mix well
2. Overlay 1 ml (approximately) of Reagent 4(Inert oil) to the above mixture
3. Close the vial immediately and incubate at 37° c.
4. Observe the change in the original blue colour of the working reagent to the brick red colour due to the hemolysate

### OBSERVATION:

1. Observe the change in colour of the reaction mixture after 30 minutes.
2. If the sample does not show decolourization, note for change in colour for every 5 -10 minutes, till the decolourization is complete.
3. If the sample does not decolourise even after 60minutes, observe the change in absorbance every 30 minutes and follow up to 4-8 hours.
4. Samples deficient in G -6-PDH may decolourise after 2-24 hours or may not decolourise at all.
5. Some samples may recolourise after decolourization , this should be ignored, the initial decolourization time should be noted
6. In case of heterozygous males or females who are carrier it is advisable to quantitatively estimate G -6-PD activity.

### REFERENCE VALUES:

Normal Range: Decolourization time (at 37° c , Hb content 15 gms/dl) : 30-60 minutes, it is recommended that each laboratory establish its own reference values.

### BIBLIOGRAPHY:

- BEUTLER, E, BLUME, K.G., KAPLAN, J.C., LOHAR, G.W., RAMOT, B. and VALENTINE, W. N.(1979)  
International Committee for Standardization in Hematology.
- DACIE V., LEWIS S. In., Practical Hematology, 7<sup>th</sup> Edition (1991), pg.204-212.

**Angstrom Biotech Pvt.Ltd.**  
G1 - 1035, RIICO Industrial Area  
Phase - III, Bhiwadi,  
Alwar, Rajasthan.  
Pin Code - 302019  
Email: info@angstrombiotech.in  
Website: www.angstrombiotech.in  
Customer Care Number - 9599194831

	Attention, see instructions for use		Consult Instructions For Use
	For in vitro diagnostic use only		Catalog #
	Store at RT		Lot Number
	Do not use if package is damaged		Date of Manufacturing
	Manufacturer		Use by